

## Product datasheet for **TR312381**

### hnRNP A1 (HNRNPA1) Human shRNA Plasmid Kit (Locus ID 3178)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	hnRNP A1 (HNRNPA1) Human shRNA Plasmid Kit (Locus ID 3178)
Locus ID:	3178
Synonyms:	ALS19; ALS20; hnRNP-A1; hnRNP A1; HNRPA1; HNRPA1L3; IBMPFD3; UP 1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	HNRNPA1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 3178). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_002136</a> , <a href="#">NM_031157</a> , <a href="#">NR_135167</a> , <a href="#">NM_002136.1</a> , <a href="#">NM_002136.3</a> , <a href="#">NM_031157.1</a> , <a href="#">NM_031157.2</a> , <a href="#">NM_031157.3</a> , <a href="#">BC009600</a> , <a href="#">BC009600.1</a> , <a href="#">BC002355</a> , <a href="#">BC004945</a> , <a href="#">BC012158</a> , <a href="#">BC020442</a> , <a href="#">BC033714</a> , <a href="#">BC035253</a> , <a href="#">BC052296</a> , <a href="#">BC070315</a> , <a href="#">BC071945</a> , <a href="#">BC073162</a> , <a href="#">BC074502</a> , <a href="#">BC078165</a> , <a href="#">BC103707</a> , <a href="#">BC104236</a> , <a href="#">BC104237</a> , <a href="#">BC121133</a> , <a href="#">NM_002136.4</a> , <a href="#">NM_031157.4</a>
UniProt ID:	<a href="#">P09651</a>
Summary:	This gene encodes a member of a family of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs), which are RNA-binding proteins that associate with pre-mRNAs in the nucleus and influence pre-mRNA processing, as well as other aspects of mRNA metabolism and transport. The protein encoded by this gene is one of the most abundant core proteins of hnRNP complexes and plays a key role in the regulation of alternative splicing. Mutations in this gene have been observed in individuals with amyotrophic lateral sclerosis 20. Multiple alternatively spliced transcript variants have been found. There are numerous pseudogenes of this gene distributed throughout the genome. [provided by RefSeq, Feb 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).