

## Product datasheet for **TR312495**

### HDAC2 Human shRNA Plasmid Kit (Locus ID 3066)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	HDAC2 Human shRNA Plasmid Kit (Locus ID 3066)
Locus ID:	3066
Synonyms:	HD2; KDAC2; RPD3; YAF1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	HDAC2 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 3066). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001527</a> , <a href="#">NR_033441</a> , <a href="#">NR_073443</a> , <a href="#">NM_001527.1</a> , <a href="#">NM_001527.2</a> , <a href="#">NM_001527.3</a> , <a href="#">BC013141</a> , <a href="#">BC031055</a> , <a href="#">BC146376</a> , <a href="#">BC148797</a> , <a href="#">BM550464</a>
UniProt ID:	<a href="#">Q92769</a>
Summary:	This gene product belongs to the histone deacetylase family. Histone deacetylases act via the formation of large multiprotein complexes, and are responsible for the deacetylation of lysine residues at the N-terminal regions of core histones (H2A, H2B, H3 and H4). This protein forms transcriptional repressor complexes by associating with many different proteins, including YY1, a mammalian zinc-finger transcription factor. Thus, it plays an important role in transcriptional regulation, cell cycle progression and developmental events. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).