

## Product datasheet for **TR312686**

### Glypican 3 (GPC3) Human shRNA Plasmid Kit (Locus ID 2719)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Glypican 3 (GPC3) Human shRNA Plasmid Kit (Locus ID 2719)
Locus ID:	2719
Synonyms:	DGSX; GTR2-2; MXR7; OCI-5; SDYS; SGB; SGBS; SGBS1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	GPC3 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 2719). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001164617</a> , <a href="#">NM_001164618</a> , <a href="#">NM_001164619</a> , <a href="#">NM_004484</a> , <a href="#">NM_004484.1</a> , <a href="#">NM_004484.2</a> , <a href="#">NM_004484.3</a> , <a href="#">NM_001164619.1</a> , <a href="#">NM_001164618.1</a> , <a href="#">NM_001164617.1</a> , <a href="#">BC035972</a> , <a href="#">BC035972.1</a> , <a href="#">NM_001164618.2</a> , <a href="#">NM_001164619.2</a> , <a href="#">NM_001164617.2</a> , <a href="#">NM_004484.4</a>
UniProt ID:	<a href="#">P51654</a>
Summary:	Cell surface heparan sulfate proteoglycans are composed of a membrane-associated protein core substituted with a variable number of heparan sulfate chains. Members of the glypican-related integral membrane proteoglycan family (GRIPS) contain a core protein anchored to the cytoplasmic membrane via a glycosyl phosphatidylinositol linkage. These proteins may play a role in the control of cell division and growth regulation. The protein encoded by this gene can bind to and inhibit the dipeptidyl peptidase activity of CD26, and it can induce apoptosis in certain cell types. Deletion mutations in this gene are associated with Simpson-Golabi-Behmel syndrome, also known as Simpson dysmorphia syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2009]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).