

Product datasheet for **TR312901**

FUSIP1 (SRSF10) Human shRNA Plasmid Kit (Locus ID 10772)

Product data:

| | |
|----------------------------------|--|
| Product Type: | shRNA Plasmids |
| Product Name: | FUSIP1 (SRSF10) Human shRNA Plasmid Kit (Locus ID 10772) |
| Locus ID: | 10772 |
| Synonyms: | FUSIP1; FUSIP2; NSSR; PPP1R149; SFRS13; SFRS13A; SRp38; SRrp40; TASR; TASR1; TASR2 |
| Vector: | pRS (TR20003) |
| E. coli Selection: | Ampicillin |
| Mammalian Cell Selection: | Puromycin |
| Format: | Retroviral plasmids |
| Components: | SRSF10 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 10772). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free. |
| RefSeq: | NM_001191005 , NM_001191006 , NM_001191007 , NM_001191009 , NM_001300936 , NM_001300937 , NM_006625 , NM_054016 , NR_034035 , NM_054016.1 , NM_054016.2 , NM_054016.3 , NM_006625.1 , NM_006625.2 , NM_006625.3 , NM_006625.4 , NM_006625.5 , NM_001191006.1 , NM_001191006.2 , NM_001191007.1 , NM_001191007.2 , NM_001191005.1 , NM_001191005.2 , NM_001191009.1 , NM_001191009.2 , NM_001300936.1 , NM_001300937.1 , BC001107 , BC010074 , BC005039 , BC071575 , NM_001300936.2 , NM_054016.4 , NM_001300937.2 , NM_001191006.3 , NM_001191005.3 , NM_001191007.3 , NM_006625.6 , NM_001191009.3 |
| UniProt ID: | O75494 |
| Summary: | This gene product is a member of the serine-arginine (SR) family of proteins, which are involved in constitutive and regulated RNA splicing. Members of this family are characterized by N-terminal RNP1 and RNP2 motifs, which are required for binding to RNA, and multiple C-terminal SR/RS repeats, which are important in mediating association with other cellular proteins. This protein interacts with the oncoprotein TLS, and abrogates the influence of TLS on adenovirus E1A pre-mRNA splicing. This gene has pseudogenes on chromosomes 4, 9, 14, 18, and 20. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014] |



[View online »](#)

shRNA Design:

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).