

## Product datasheet for **TR312951**

### FNTA Human shRNA Plasmid Kit (Locus ID 2339)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	FNTA Human shRNA Plasmid Kit (Locus ID 2339)
Locus ID:	2339
Synonyms:	FPTA; MGC99680; PGGT1A; PTAR2
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	FNTA - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 2339). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001018676</a> , <a href="#">NM_001018677</a> , <a href="#">NM_002027</a> , <a href="#">NR_033698</a> , <a href="#">NM_002027.1</a> , <a href="#">NM_002027.2</a> , <a href="#">NM_001018676.1</a> , <a href="#">NM_001018677.1</a> , <a href="#">BC084566</a> , <a href="#">BC084566.1</a> , <a href="#">BC017029</a> , <a href="#">BC037295</a> , <a href="#">BM686906</a> , <a href="#">NM_002027.3</a>
UniProt ID:	<a href="#">P49354</a>
Summary:	Prenyltransferases can attach either a farnesyl group or a geranylgeranyl group in thioether linkage to the cysteine residue of proteins with a C-terminal CAAX box. CAAX geranylgeranyltransferase and CAAX farnesyltransferase are heterodimers that share the same alpha subunit but have different beta subunits. This gene encodes the alpha subunit of these transferases. Alternative splicing results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 11 and 13. [provided by RefSeq, May 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).