

## Product datasheet for **TR313152**

### ER81 (ETV1) Human shRNA Plasmid Kit (Locus ID 2115)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	ER81 (ETV1) Human shRNA Plasmid Kit (Locus ID 2115)
Locus ID:	2115
Synonyms:	ER81
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	ETV1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 2115). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001163147</a> , <a href="#">NM_001163148</a> , <a href="#">NM_001163149</a> , <a href="#">NM_001163150</a> , <a href="#">NM_001163151</a> , <a href="#">NM_001163152</a> , <a href="#">NM_004956</a> , <a href="#">NR_120445</a> , <a href="#">NM_004956.1</a> , <a href="#">NM_004956.2</a> , <a href="#">NM_004956.3</a> , <a href="#">NM_004956.4</a> , <a href="#">NM_001163152.1</a> , <a href="#">NM_001163151.1</a> , <a href="#">NM_001163150.1</a> , <a href="#">NM_001163147.1</a> , <a href="#">NM_001163148.1</a> , <a href="#">NM_001163149.1</a> , <a href="#">BC098403</a> , <a href="#">BC098403.1</a> , <a href="#">BC042950</a> , <a href="#">BC045776</a> , <a href="#">BC106762</a> , <a href="#">BC106763</a> , <a href="#">NM_001370556</a> , <a href="#">NM_001370555</a> , <a href="#">NM_004956.5</a>
UniProt ID:	<a href="#">P50549</a>
Summary:	This gene encodes a member of the ETS (E twenty-six) family of transcription factors. The ETS proteins regulate many target genes that modulate biological processes like cell growth, angiogenesis, migration, proliferation and differentiation. All ETS proteins contain an ETS DNA-binding domain that binds to DNA sequences containing the consensus 5'-CGGA[AT]-3'. The protein encoded by this gene contains a conserved short acidic transactivation domain (TAD) in the N-terminal region, in addition to the ETS DNA-binding domain in the C-terminal region. This gene is involved in chromosomal translocations, which result in multiple fusion proteins including EWS-ETV1 in Ewing sarcoma and at least 10 ETV1 partners (see PMID: 19657377, Table 1) in prostate cancer. In addition to chromosomal rearrangement, this gene is overexpressed in prostate cancer, melanoma and gastrointestinal stromal tumor. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2016]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).