

Product datasheet for **TR313406**

Dynamin 2 (DNM2) Human shRNA Plasmid Kit (Locus ID 1785)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Dynamin 2 (DNM2) Human shRNA Plasmid Kit (Locus ID 1785)
Locus ID:	1785
Synonyms:	CMT2M; CMTDI1; CMTDIB; DI-CMTB; DYN2; DYNII; LCCS5
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	DNM2 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 1785). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001005360 , NM_001005361 , NM_001005362 , NM_001190716 , NM_004945 , NM_004945.1 , NM_004945.2 , NM_004945.3 , NM_001005360.1 , NM_001005360.2 , NM_001005361.1 , NM_001005361.2 , NM_001005362.1 , NM_001005362.2 , NM_001190716.1 , BC054501 , BC054501.1 , BC039596 , BM148999 , NM_001005361.3 , NM_001190716.2 , NM_001005360.3
UniProt ID:	P50570
Summary:	Dynamins represent one of the subfamilies of GTP-binding proteins. These proteins share considerable sequence similarity over the N-terminal portion of the molecule, which contains the GTPase domain. Dynamins are associated with microtubules. They have been implicated in cell processes such as endocytosis and cell motility, and in alterations of the membrane that accompany certain activities such as bone resorption by osteoclasts. Dynamins bind many proteins that bind actin and other cytoskeletal proteins. Dynamins can also self-assemble, a process that stimulates GTPase activity. Five alternatively spliced transcripts encoding different proteins have been described. Additional alternatively spliced transcripts may exist, but their full-length nature has not been determined. [provided by RefSeq, Jun 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).