

## Product datasheet for **TR313478**

### DHX30 Human shRNA Plasmid Kit (Locus ID 22907)

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | shRNA Plasmids   |
| Product Name:             | DHX30 Human shRNA Plasmid Kit (Locus ID 22907)   |
| Locus ID:                 | 22907  |
| Synonyms:                 | DDX30; NEDMIAL; RETCOR   |
| Vector:                   | pRS (TR20003)  |
| E. coli Selection:        | Ampicillin   |
| Mammalian Cell Selection: | Puromycin  |
| Format:                   | Retroviral plasmids  |
| Components:               | DHX30 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 22907). 5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.  |
| RefSeq:                   | <a href="#">BC015029</a> , <a href="#">NM_001330990</a> , <a href="#">NM_014966</a> , <a href="#">NM_138614</a> , <a href="#">NM_138615</a> , <a href="#">NR_075079</a> , <a href="#">NM_138615.1</a> , <a href="#">NM_138615.2</a> , <a href="#">NM_014966.1</a> , <a href="#">NM_014966.2</a> , <a href="#">NM_014966.3</a> , <a href="#">BC015029.1</a> , <a href="#">NM_138614.1</a> , <a href="#">BC020126</a> , <a href="#">BC020126.1</a> , <a href="#">BC038417</a> , <a href="#">BC038417.1</a> , <a href="#">BC014237</a> , <a href="#">BC028405</a> , <a href="#">BC047335</a> , <a href="#">BC063510</a> , <a href="#">NM_138615.3</a>   |
| UniProt ID:               | <a href="#">Q7L2E3</a>   |
| Summary:                  | DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of this DEAD box protein family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. The family member encoded by this gene is a mitochondrial nucleoid protein that associates with mitochondrial DNA. It has also been identified as a component of a transcriptional repressor complex that functions in retinal development, and it is required to optimize the function of the zinc-finger antiviral protein. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Feb 2013] |



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).