

## Product datasheet for **TR313875**

### CLCNKA Human shRNA Plasmid Kit (Locus ID 1187)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CLCNKA Human shRNA Plasmid Kit (Locus ID 1187)
Locus ID:	1187
Synonyms:	CLC-K1; CLCK1; hCLC-Ka
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	CLCNKA - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 1187). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001042704</a> , <a href="#">NM_001257139</a> , <a href="#">NM_004070</a> , <a href="#">NM_004070.1</a> , <a href="#">NM_004070.2</a> , <a href="#">NM_004070.3</a> , <a href="#">NM_001042704.1</a> , <a href="#">NM_001257139.1</a> , <a href="#">BC053869</a> , <a href="#">BC053869.1</a> , <a href="#">BC035373</a> , <a href="#">BC048282</a> , <a href="#">NM_004070.4</a> , <a href="#">NM_001042704.2</a> , <a href="#">NM_001257139.2</a>
UniProt ID:	<a href="#">P51800</a>
Summary:	This gene is a member of the CLC family of voltage-gated chloride channels. The encoded protein is predicted to have 12 transmembrane domains, and requires a beta subunit called barttin to form a functional channel. It is thought to function in salt reabsorption in the kidney and potassium recycling in the inner ear. The gene is highly similar to CLCNKB, which is located 10 kb downstream from this gene. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).