

## Product datasheet for **TR313977**

### CEP57 Human shRNA Plasmid Kit (Locus ID 9702)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CEP57 Human shRNA Plasmid Kit (Locus ID 9702)
Locus ID:	9702
Synonyms:	MVA2; PIG8; TSP57
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	CEP57 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 9702). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC009053</a> , <a href="#">NM_001243776</a> , <a href="#">NM_001243777</a> , <a href="#">NM_014679</a> , <a href="#">NM_014679.1</a> , <a href="#">NM_014679.2</a> , <a href="#">NM_014679.3</a> , <a href="#">NM_014679.4</a> , <a href="#">NM_001243777.1</a> , <a href="#">NM_001243776.1</a> , <a href="#">BC001233</a> , <a href="#">BC029385</a> , <a href="#">BC031048</a> , <a href="#">BC039711</a> , <a href="#">NM_001363604</a> , <a href="#">NM_014679.5</a>
UniProt ID:	<a href="#">Q86XR8</a>
Summary:	This gene encodes a cytoplasmic protein called Translokin. This protein localizes to the centrosome and has a function in microtubular stabilization. The N-terminal half of this protein is required for its centrosome localization and for its multimerization, and the C-terminal half is required for nucleating, bundling and anchoring microtubules to the centrosomes. This protein specifically interacts with fibroblast growth factor 2 (FGF2), sorting nexin 6, Ran-binding protein M and the kinesins KIF3A and KIF3B, and thus mediates the nuclear translocation and mitogenic activity of the FGF2. It also interacts with cyclin D1 and controls nucleocytoplasmic distribution of the cyclin D1 in quiescent cells. This protein is crucial for maintaining correct chromosomal number during cell division. Mutations in this gene cause mosaic variegated aneuploidy syndrome, a rare autosomal recessive disorder. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Aug 2011]



[View online »](#)

**shRNA Design:**

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

**Performance Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).