

## Product datasheet for **TR314196**

### CAPZB Human shRNA Plasmid Kit (Locus ID 832)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CAPZB Human shRNA Plasmid Kit (Locus ID 832)
Locus ID:	832
Synonyms:	CAPB; CAPPB; CAPZ
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	CAPZB - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 832). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001206540</a> , <a href="#">NM_001206541</a> , <a href="#">NM_001282162</a> , <a href="#">NM_001313932</a> , <a href="#">NM_004930</a> , <a href="#">NR_038125</a> , <a href="#">NM_004930.1</a> , <a href="#">NM_004930.2</a> , <a href="#">NM_004930.3</a> , <a href="#">NM_004930.4</a> , <a href="#">NM_001206540.1</a> , <a href="#">NM_001206540.2</a> , <a href="#">NM_001206541.1</a> , <a href="#">NM_001206541.2</a> , <a href="#">NM_001282162.1</a> , <a href="#">BC107752</a> , <a href="#">BC107752.1</a> , <a href="#">BC008095</a> , <a href="#">BC012305</a> , <a href="#">BC024601</a> , <a href="#">BC069031</a> , <a href="#">BC109241</a> , <a href="#">BC109242</a> , <a href="#">NM_001206541.3</a> , <a href="#">NM_001206540.3</a> , <a href="#">NM_004930.5</a> , <a href="#">NM_001282162.2</a>
UniProt ID:	<a href="#">P47756</a>
Summary:	This gene encodes the beta subunit of the barbed-end actin binding protein, which belongs to the F-actin capping protein family. The capping protein is a heterodimeric actin capping protein that blocks actin filament assembly and disassembly at the fast growing (barbed) filament ends and functions in regulating actin filament dynamics as well as in stabilizing actin filament lengths in muscle and nonmuscle cells. A pseudogene of this gene is located on the long arm of chromosome 2. Multiple alternatively spliced transcript variants encoding different isoforms have been found.[provided by RefSeq, Aug 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).