

Product datasheet for **TR314529**

Beta 1,4 galactosyltransferase 6 (B4GALT6) Human shRNA Plasmid Kit (Locus ID 9331)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Beta 1,4 galactosyltransferase 6 (B4GALT6) Human shRNA Plasmid Kit (Locus ID 9331)
Locus ID:	9331
Synonyms:	B4Gal-T6; beta4Gal-T6
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	B4GALT6 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 9331). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_004775 , NM_001330570 , NM_004775.1 , NM_004775.2 , NM_004775.3 , BC069642 , BC069642.1 , BC069620 , BC074835 , BC074884 , NM_004775.5
UniProt ID:	Q9UBX8
Summary:	This gene is one of seven beta-1,4-galactosyltransferase (beta4GalT) genes in human. They encode type II membrane-bound glycoproteins that appear to have exclusive specificity for the donor substrate UDP-galactose; all transfer galactose in a beta1,4 linkage to similar acceptor sugars: GlcNAc, Glc, and Xyl. Each beta4GalT has a distinct function in the biosynthesis of different glycoconjugates and saccharide structures. As type II membrane proteins, they have an N-terminal hydrophobic signal sequence that directs the protein to the Golgi apparatus and which then remains uncleaved to function as a transmembrane anchor. This gene produces multiple protein isoforms - some of which are predicted to lack the N-terminal hydrophobic signal sequence and transmembrane domain. By sequence similarity, the beta4GalTs form four groups: beta4GalT1 and beta4GalT2, beta4GalT3 and beta4GalT4, beta4GalT5 and beta4GalT6, and beta4GalT7. The canonical enzyme encoded by this gene is a lactosylceramide synthase important for glycolipid biosynthesis. [provided by RefSeq, Jan 2020]



[View online »](#)

shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .
Performance Guaranteed:	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>