

## Product datasheet for **TR314702**

### ARF1 Human shRNA Plasmid Kit (Locus ID 375)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	ARF1 Human shRNA Plasmid Kit (Locus ID 375)
Locus ID:	375
Synonyms:	PVNH8
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	ARF1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 375). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001024226</a> , <a href="#">NM_001024227</a> , <a href="#">NM_001024228</a> , <a href="#">NM_001658</a> , <a href="#">NM_001024227.1</a> , <a href="#">NM_001658.1</a> , <a href="#">NM_001658.2</a> , <a href="#">NM_001658.3</a> , <a href="#">NM_001024228.1</a> , <a href="#">NM_001024226.1</a> , <a href="#">BC009247</a> , <a href="#">BC009247.2</a> , <a href="#">BC011358</a> , <a href="#">BC011358.1</a> , <a href="#">BC010415</a> , <a href="#">BC010429</a> , <a href="#">BM353046</a> , <a href="#">BM844480</a> , <a href="#">NM_001024228.2</a> , <a href="#">NM_001658.4</a>
UniProt ID:	<a href="#">P84077</a>
Summary:	ADP-ribosylation factor 1 (ARF1) is a member of the human ARF gene family. The family members encode small guanine nucleotide-binding proteins that stimulate the ADP-ribosyltransferase activity of cholera toxin and play a role in vesicular trafficking as activators of phospholipase D. The gene products, including 6 ARF proteins and 11 ARF-like proteins, constitute a family of the RAS superfamily. The ARF proteins are categorized as class I (ARF1, ARF2 and ARF3), class II (ARF4 and ARF5) and class III (ARF6), and members of each class share a common gene organization. The ARF1 protein is localized to the Golgi apparatus and has a central role in intra-Golgi transport. Multiple alternatively spliced transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).