

## Product datasheet for **TR314717**

### APOL3 Human shRNA Plasmid Kit (Locus ID 80833)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	APOL3 Human shRNA Plasmid Kit (Locus ID 80833)
Locus ID:	80833
Synonyms:	apoL-III; APOLIII; CG12_1; CG121
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	APOL3 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 80833). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_014349</a> , <a href="#">NM_030644</a> , <a href="#">NM_145639</a> , <a href="#">NM_145640</a> , <a href="#">NM_145641</a> , <a href="#">NM_145642</a> , <a href="#">NR_027833</a> , <a href="#">NR_027834</a> , <a href="#">NR_027835</a> , <a href="#">NM_014349.1</a> , <a href="#">NM_014349.2</a> , <a href="#">NM_145639.1</a> , <a href="#">NM_145641.1</a> , <a href="#">NM_145641.2</a> , <a href="#">NM_145642.1</a> , <a href="#">NM_145642.2</a> , <a href="#">NM_030644.1</a> , <a href="#">NM_145640.1</a> , <a href="#">NM_145640.2</a> , <a href="#">BC042918</a> , <a href="#">BC042918.1</a> , <a href="#">BC050596</a> , <a href="#">BC038352</a> , <a href="#">BC070144</a> , <a href="#">NM_030644.2</a> , <a href="#">NM_145639.2</a> , <a href="#">NM_014349.3</a>
UniProt ID:	<a href="#">O95236</a>
Summary:	This gene is a member of the apolipoprotein L gene family, and it is present in a cluster with other family members on chromosome 22. The encoded protein is found in the cytoplasm, where it may affect the movement of lipids, including cholesterol, and/or allow the binding of lipids to organelles. In addition, expression of this gene is up-regulated by tumor necrosis factor-alpha in endothelial cells lining the normal and atherosclerotic iliac artery and aorta. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).