

Product datasheet for TR314737

APC Human shRNA Plasmid Kit (Locus ID 324)

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	shRNA Plasmids
Product Name:	APC Human shRNA Plasmid Kit (Locus ID 324)
Locus ID:	324
Synonyms:	BTPS2; DESMD; DP2; DP2.5; DP3; GS; PPP1R46
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	APC - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 324). 5μg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<u>NM 000038, NM 001127510, NM 001127511, NM 001354895, NM 001354896,</u> <u>NM 001354897, NM 001354898, NM 001354899, NM 001354900, NM 001354901,</u> <u>NM 001354902, NM 001354903, NM 001354904, NM 001354905, NM 001354906,</u> <u>NM 000038.1, NM 000038.2, NM 000038.3, NM 000038.4, NM 000038.5, NM 001127510.1,</u> <u>NM 001127510.2, NM 001127511.1, NM 001127511.2, BC034955, BC056268, BC111462,</u> <u>BC111591, BC111930, BM271918, BM713901, BM975470, NM 000038.6, NM 001127510.3</u>
UniProt ID:	<u>P25054</u>
Summary:	This gene encodes a tumor suppressor protein that acts as an antagonist of the Wnt signaling pathway. It is also involved in other processes including cell migration and adhesion, transcriptional activation, and apoptosis. Defects in this gene cause familial adenomatous polyposis (FAP), an autosomal dominant pre-malignant disease that usually progresses to malignancy. Mutations in the APC gene have been found to occur in most colorectal cancers. Disease-associated mutations tend to be clustered in a small region designated the mutation cluster region (MCR) and result in a truncated protein product. [provided by RefSeq, Dec 2019]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u> . If you need a special design or shRNA sequence, please utilize our <u>custom shRNA service</u> .



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

GRIGENE APC Human shRNA Plasmid Kit (Locus ID 324) – TR314737

Performance Guaranteed: OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., <u>9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US</u>