

## Product datasheet for **TR314797**

### ANGPTL1 Human shRNA Plasmid Kit (Locus ID 9068)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	ANGPTL1 Human shRNA Plasmid Kit (Locus ID 9068)
Locus ID:	9068
Synonyms:	ANG3; ANGPT3; AngY; ARP1; dj595C2.2; UNQ162
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	ANGPTL1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 9068). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_004673</a> , <a href="#">NM_004673.1</a> , <a href="#">NM_004673.2</a> , <a href="#">NM_004673.3</a> , <a href="#">BC050640</a> , <a href="#">BC050640.1</a> , <a href="#">BC047670</a> , <a href="#">BM710120</a> , <a href="#">NM_004673.4</a>
UniProt ID:	<a href="#">O95841</a>
Summary:	Angiopoietins are members of the vascular endothelial growth factor family and the only known growth factors largely specific for vascular endothelium. Angiopoietin-1, angiopoietin-2, and angiopoietin-4 participate in the formation of blood vessels. The protein encoded by this gene is another member of the angiopoietin family that is widely expressed in adult tissues with mRNA levels highest in highly vascularized tissues. This protein was found to be a secretory protein that does not act as an endothelial cell mitogen in vitro. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).