

Product datasheet for **TR315009**

ABR Human shRNA Plasmid Kit (Locus ID 29)

Product data:

Product Type:	shRNA Plasmids
Product Name:	ABR Human shRNA Plasmid Kit (Locus ID 29)
Locus ID:	29
Synonyms:	MDB
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	ABR - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 29). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001092 , NM_001159746 , NM_001256847 , NM_001282149 , NM_001322840 , NM_001322841 , NM_001322842 , NM_021962 , NR_046393 , NM_001092.1 , NM_001092.2 , NM_001092.3 , NM_001092.4 , NM_021962.1 , NM_021962.2 , NM_021962.3 , NM_021962.4 , NM_001159746.1 , NM_001159746.2 , NM_001256847.1 , NM_001256847.2 , NM_001282149.1 , BC036807 , BC141909 , BC156483 , BC172517 , BM547076 , BM827672 , NM_001092.5 , NM_001282149.2 , NM_001159746.3 , NM_001256847.3 , NM_021962.5
UniProt ID:	Q12979
Summary:	This gene encodes a protein that is similar to the protein encoded by the breakpoint cluster region gene located on chromosome 22. The protein encoded by this gene contains a GTPase-activating protein domain, a domain found in members of the Rho family of GTP-binding proteins. Functional studies in mice determined that this protein plays a role in vestibular morphogenesis. Alternatively spliced transcript variants have been reported for this gene. [provided by RefSeq, Feb 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).