

## Product datasheet for **TR315487**

### IRF7 Human shRNA Plasmid Kit (Locus ID 3665)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	IRF7 Human shRNA Plasmid Kit (Locus ID 3665)
Locus ID:	3665
Synonyms:	IMD39; IRF-7; IRF-7H; IRF7A; IRF7B; IRF7C; IRF7H
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	IRF7 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 3665). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001572</a> , <a href="#">NM_004029</a> , <a href="#">NM_004030</a> , <a href="#">NM_004031</a> , <a href="#">NM_004029.1</a> , <a href="#">NM_004029.2</a> , <a href="#">NM_001572.1</a> , <a href="#">NM_001572.2</a> , <a href="#">NM_001572.3</a> , <a href="#">NM_004031.1</a> , <a href="#">BC136555</a> , <a href="#">BC144242</a> , <a href="#">BM849172</a> , <a href="#">NM_004029.4</a> , <a href="#">NM_001572.5</a>
UniProt ID:	<a href="#">Q92985</a>
Summary:	IRF7 encodes interferon regulatory factor 7, a member of the interferon regulatory transcription factor (IRF) family. IRF7 has been shown to play a role in the transcriptional activation of virus-inducible cellular genes, including interferon beta chain genes. Inducible expression of IRF7 is largely restricted to lymphoid tissue. Multiple IRF7 transcript variants have been identified, although the functional consequences of these have not yet been established. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).