

## Product datasheet for **TR316708**

### RPS14 Human shRNA Plasmid Kit (Locus ID 6208)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	RPS14 Human shRNA Plasmid Kit (Locus ID 6208)
Locus ID:	6208
Synonyms:	EMTB; S14
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	RPS14 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 6208). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001025070</a> , <a href="#">NM_001025071</a> , <a href="#">NM_005617</a> , <a href="#">NM_001025070.1</a> , <a href="#">NM_005617.1</a> , <a href="#">NM_005617.2</a> , <a href="#">NM_005617.3</a> , <a href="#">NM_001025071.1</a> , <a href="#">BC020515</a> , <a href="#">BC020515.1</a> , <a href="#">BC091474</a> , <a href="#">BC091474.1</a> , <a href="#">BC001126</a> , <a href="#">BC003401</a> , <a href="#">BC006784</a> , <a href="#">BC110792</a> , <a href="#">BM454896</a> , <a href="#">BM557014</a> , <a href="#">NM_001025071.2</a>
UniProt ID:	<a href="#">P62263</a>
Summary:	Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S11P family of ribosomal proteins. It is located in the cytoplasm. Transcript variants utilizing alternative transcription initiation sites have been described in the literature. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. In Chinese hamster ovary cells, mutations in this gene can lead to resistance to emetine, a protein synthesis inhibitor. Multiple alternatively spliced transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).