

## Product datasheet for **TR316769**

### DAP12 (TYROBP) Human shRNA Plasmid Kit (Locus ID 7305)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	DAP12 (TYROBP) Human shRNA Plasmid Kit (Locus ID 7305)
Locus ID:	7305
Synonyms:	DAP12; KARAP; PLOSL; PLOSL1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	TYROBP - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 7305). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001173514</a> , <a href="#">NM_001173515</a> , <a href="#">NM_003332</a> , <a href="#">NM_198125</a> , <a href="#">NR_033390</a> , <a href="#">NM_198125.1</a> , <a href="#">NM_198125.2</a> , <a href="#">NM_003332.1</a> , <a href="#">NM_003332.2</a> , <a href="#">NM_003332.3</a> , <a href="#">NM_001173515.1</a> , <a href="#">NM_001173514.1</a> , <a href="#">BC011175</a> , <a href="#">BC011175.1</a> , <a href="#">NM_003332.4</a> , <a href="#">NM_001173515.2</a> , <a href="#">NM_001173514.2</a>
UniProt ID:	<a href="#">O43914</a>
Summary:	This gene encodes a transmembrane signaling polypeptide which contains an immunoreceptor tyrosine-based activation motif (ITAM) in its cytoplasmic domain. The encoded protein may associate with the killer-cell inhibitory receptor (KIR) family of membrane glycoproteins and may act as an activating signal transduction element. This protein may bind zeta-chain (TCR) associated protein kinase 70kDa (ZAP-70) and spleen tyrosine kinase (SYK) and play a role in signal transduction, bone modeling, brain myelination, and inflammation. Mutations within this gene have been associated with polycystic lipomembranous osteodysplasia with sclerosing leukoencephalopathy (PLOSL), also known as Nasu-Hakola disease. Its putative receptor, triggering receptor expressed on myeloid cells 2 (TREM2), also causes PLOSL. Multiple alternative transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Mar 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).