

Product datasheet for **TR316997**

alpha Glucosidase II (GANAB) Human shRNA Plasmid Kit (Locus ID 23193)

Product data:

Product Type:	shRNA Plasmids
Product Name:	alpha Glucosidase II (GANAB) Human shRNA Plasmid Kit (Locus ID 23193)
Locus ID:	23193
Synonyms:	G2AN; GIIA; GLUII; PKD3
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	GANAB - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 23193). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001278192 , NM_001278193 , NM_001278194 , NM_001329222 , NM_001329223 , NM_001329224 , NM_001329225 , NM_014610 , NM_198334 , NM_198335 , NM_198334.1 , NM_198334.2 , NM_198335.1 , NM_198335.2 , NM_198335.3 , NM_001278193.1 , NM_001278194.1 , NM_001278192.1 , BC005405 , BC017433 , BC017435 , BC034439 , BC065266 , NM_001278193.2 , NM_198334.3 , NM_198335.4 , NM_001278192.2
UniProt ID:	Q14697
Summary:	This gene encodes the alpha subunit of glucosidase II and a member of the glycosyl hydrolase 31 family of proteins. The heterodimeric enzyme glucosidase II plays a role in protein folding and quality control by cleaving glucose residues from immature glycoproteins in the endoplasmic reticulum. Expression of the encoded protein is elevated in lung tumor tissue and in response to UV irradiation. Mutations in this gene cause autosomal-dominant polycystic kidney and liver disease. [provided by RefSeq, Jul 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).