

## Product datasheet for **TR317021**

### HP1 alpha (CBX5) Human shRNA Plasmid Kit (Locus ID 23468)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	HP1 alpha (CBX5) Human shRNA Plasmid Kit (Locus ID 23468)
Locus ID:	23468
Synonyms:	HEL25; HP1; HP1A
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	CBX5 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 23468). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001127321</a> , <a href="#">NM_001127322</a> , <a href="#">NM_012117</a> , <a href="#">NM_012117.1</a> , <a href="#">NM_012117.2</a> , <a href="#">NM_001127321.1</a> , <a href="#">NM_001127322.1</a> , <a href="#">BC006821</a> , <a href="#">BC006821.1</a> , <a href="#">BC030149</a> , <a href="#">BC047432</a> , <a href="#">NM_012117.3</a>
UniProt ID:	<a href="#">P45973</a>
Summary:	This gene encodes a highly conserved nonhistone protein, which is a member of the heterochromatin protein family. The protein is enriched in the heterochromatin and associated with centromeres. The protein has a single N-terminal chromodomain which can bind to histone proteins via methylated lysine residues, and a C-terminal chromo shadow-domain (CSD) which is responsible for the homodimerization and interaction with a number of chromatin-associated nonhistone proteins. The encoded product is involved in the formation of functional kinetochore through interaction with essential kinetochore proteins. The gene has a pseudogene located on chromosome 3. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).