

## Product datasheet for **TR317442**

### LYNX1 Human shRNA Plasmid Kit (Locus ID 66004)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	LYNX1 Human shRNA Plasmid Kit (Locus ID 66004)
Locus ID:	66004
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	LYNX1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 66004). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_177457</a> , <a href="#">NM_177476</a> , <a href="#">NM_177477</a> , <a href="#">NM_001356370</a> , <a href="#">NM_023946.1</a> , <a href="#">NM_023946.2</a> , <a href="#">NM_023946.3</a> , <a href="#">NM_177458.1</a> , <a href="#">NM_177476.1</a> , <a href="#">NM_177476.3</a> , <a href="#">NM_177477.1</a> , <a href="#">NM_177477.2</a> , <a href="#">NM_177477.3</a> , <a href="#">NM_177457.1</a> , <a href="#">NM_177457.2</a> , <a href="#">NM_177457.3</a> , <a href="#">NM_177457.4</a> , <a href="#">BC014560</a> , <a href="#">BC047923</a> , <a href="#">BM172253</a> , <a href="#">BM671222</a> , <a href="#">NM_177457.5</a> , <a href="#">NM_177477.4</a>
UniProt ID:	<a href="#">P0DP58</a>
Summary:	This gene encodes a GPI-anchored, cell membrane bound member of the Ly6/uPAR (LU) superfamily of proteins containing the unique three-finger LU domain. This protein interacts with nicotinic acetylcholine receptors (nAChRs), and is thought to function as a modulator of nAChR activity to prevent excessive excitation. Alternatively spliced transcript variants have been found for this gene. Read-through transcription between this gene and the neighboring downstream gene (SLURP2) generates naturally-occurring transcripts (LYNX1-SLURP2) that encode a fusion protein comprised of sequence sharing identity with each individual gene product. [provided by RefSeq, Sep 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).