

Product datasheet for **TR317790**

RAB3IP Human shRNA Plasmid Kit (Locus ID 117177)

Product data:

Product Type:	shRNA Plasmids
Product Name:	RAB3IP Human shRNA Plasmid Kit (Locus ID 117177)
Locus ID:	117177
Synonyms:	RABIN3; RABIN8
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	RAB3IP - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 117177). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001024647 , NM_001278402 , NM_022456 , NM_175623 , NM_175624 , NM_175625 , NM_175626 , NM_175627 , NR_103519 , NR_103520 , NM_001024647.1 , NM_001024647.2 , NM_001024647.3 , NM_175624.1 , NM_175624.2 , NM_175624.3 , NM_022456.1 , NM_022456.2 , NM_022456.3 , NM_022456.4 , NM_175623.1 , NM_175623.2 , NM_175623.3 , NM_175625.1 , NM_175625.2 , NM_175625.3 , NM_001278402.1 , BC015548 , BC015548.1 , BC014553 , BC059358 , NM_022456.5 , NM_175623.4 , NM_175625.4 , NM_175624.4
UniProt ID:	Q96QF0
Summary:	Guanine nucleotide exchange factor (GEF) which may activate RAB8A and RAB8B. Promotes the exchange of GDP to GTP, converting inactive GDP-bound Rab proteins into their active GTP-bound form. Mediates the release of GDP from RAB8A and RAB8B but not from RAB3A or RAB5. Modulates actin organization and promotes polarized transport of RAB8A-specific vesicles to the cell surface. Together with RAB11A, RAB8A, the exocyst complex, PARD3, PRKCI, ANXA2, CDC42 and DNMBP promotes transcytosis of PODXL to the apical membrane initiation sites (AMIS), apical surface formation and lumenogenesis.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).