

Product datasheet for **TR319767**

DLEU2 Human shRNA Plasmid Kit (Locus ID 8847)

Product data:

| | |
|---------------------------|---|
| Product Type: | shRNA Plasmids |
| Product Name: | DLEU2 Human shRNA Plasmid Kit (Locus ID 8847) |
| Locus ID: | 8847 |
| Synonyms: | 1810047A16Rik; 4932411A06Rik; 7730401J12Rik; AI197277; AI788824; Alt1; deleted in lymphocytic leukemia, 2; DLB2, LEU2, BCMSUN, MGC88287; Leu2 |
| Vector: | pRS (TR20003) |
| E. coli Selection: | Ampicillin |
| Mammalian Cell Selection: | Puromycin |
| Format: | Retroviral plasmids |
| Components: | DLEU2 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 8847). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free. |
| RefSeq: | NM_006021 , NR_002612 , NR_152566 , NR_152567 , NR_152568 , NR_152569 , NR_152570 , NR_152571 , NR_152572 , NM_006021.1 , BC006995 , BC017819 , BC022282 , BC022349 , BC026275 , BC030971 , BC033120 , BC038090 , BC055417 , BC070292 |
| Summary: | This locus represents a microRNA host gene and also produces long alternatively spliced non-coding RNAs. This genome region was observed to be deleted or epigenetically suppressed in leukemia, and was implicated as a negative regulator of cell proliferation. However, an alternative transcript produced at this locus was also found to promote progression through the cell cycle via angiotensin I converting enzyme 2 and cyclin D1. [provided by RefSeq, Dec 2017] |
| shRNA Design: | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service . |



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).