

## Product datasheet for **TR320076**

### BEAN1 Human shRNA Plasmid Kit (Locus ID 146227)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	BEAN1 Human shRNA Plasmid Kit (Locus ID 146227)
Locus ID:	146227
Synonyms:	BEAN; SCA31
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	BEAN1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector (Gene ID = 146227). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001136106</a> , <a href="#">NM_001178020</a> , <a href="#">NM_001197224</a> , <a href="#">NM_001197225</a> , <a href="#">NM_001136106.1</a> , <a href="#">NM_001136106.2</a> , <a href="#">NM_001136106.3</a> , <a href="#">NM_001136106.4</a> , <a href="#">NM_001178020.1</a> , <a href="#">NM_001178020.2</a> , <a href="#">NM_001197224.1</a> , <a href="#">NM_001197224.2</a> , <a href="#">NM_001197224.3</a> , <a href="#">NM_001197225.1</a> , <a href="#">NM_001197225.2</a> , <a href="#">NM_001197225.3</a> , <a href="#">BC000818</a> , <a href="#">BC105723</a> , <a href="#">BC105724</a> , <a href="#">BC107477</a> , <a href="#">BM666201</a> , <a href="#">NM_001178020.3</a> , <a href="#">NM_001197225.4</a> , <a href="#">NM_001197224.4</a> , <a href="#">NM_001136106.5</a>
UniProt ID:	<a href="#">Q3B7T3</a>
Summary:	The protein encoded by this gene is one of several proteins that interact with NEDD4, a member of a family of ubiquitin-protein ligases. These proteins have PY motifs in common that bind to the WW domains of NEDD4. NEDD4 is developmentally regulated, and is highly expressed in embryonic tissues. Mutations in this gene (i.e., intronic insertions of >100 copies of pentanucleotide repeats including a (TGGAA) <sub>n</sub> sequence) are associated with spinocerebellar ataxia type 31. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).