

Product datasheet for **TR320203**

CCDC88C Human shRNA Plasmid Kit (Locus ID 440193)

Product data:

Product Type:	shRNA Plasmids
Product Name:	CCDC88C Human shRNA Plasmid Kit (Locus ID 440193)
Locus ID:	440193
Synonyms:	DAPLE; HKRP2; HYC1; KIAA1509; SCA40
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	CCDC88C - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 440193). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC028565 , NM_001080414 , NM_001080414.1 , NM_001080414.2 , NM_001080414.3 , BC032316 , BC035914
UniProt ID:	Q9P219
Summary:	This gene encodes a ubiquitously expressed coiled-coil domain-containing protein that interacts with the dishevelled protein and is a negative regulator of the Wnt signalling pathway. The protein encoded by this gene has a PDZ-domain binding motif in its C-terminus with which it interacts with the dishevelled protein. Dishevelled is a scaffold protein involved in the regulation of the Wnt signaling pathway. The Wnt signaling pathway plays an important role in embryonic development, tissue maintenance, and cancer progression. Mutations in this gene cause autosomal recessive, primary non-syndromic congenital hydrocephalus; a condition characterized by excessive accumulation of cerebrospinal fluid in the ventricles of the brain. [provided by RefSeq, Jan 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).