

## **Product datasheet for TR320317**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## CKII alpha (CSNK2A1) Human shRNA Plasmid Kit (Locus ID 1457)

**Product data:** 

**Product Type:** shRNA Plasmids

**Product Name:** CKII alpha (CSNK2A1) Human shRNA Plasmid Kit (Locus ID 1457)

Locus ID: 1457

Synonyms: CK2A1; Cka1; Cka2; CKII; OCNDS

**Vector:** pRS (TR20003)

E. coli Selection: Ampicillin

Mammalian Cell Puromycin

Selection: Format:

Retroviral plasmids

Components: CSNK2A1 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID =

1457). 5µg purified plasmid DNA per construct

29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.

RefSeq: NM 001895, NM 177559, NM 177560, NM 177560.1, NM 177560.2, NM 177559.1,

NM 177559.2, NM 001895.1, NM 001895.2, NM 001895.3, BC053532, BC053532.1, BC050036,

BC002615, BC011668, BC071167, NM 001362771, NM 001362770, NM 001895.4,

NM 177559.3

UniProt ID: P68400

**Summary:** Casein kinase II is a serine/threonine protein kinase that phosphorylates acidic proteins such

as casein. It is involved in various cellular processes, including cell cycle control, apoptosis, and circadian rhythm. The kinase exists as a tetramer and is composed of an alpha, an alphaprime, and two beta subunits. The alpha subunits contain the catalytic activity while the beta subunits undergo autophosphorylation. The protein encoded by this gene represents the alpha subunit. Multiple transcript variants encoding different protein isoforms have been

found for this gene. [provided by RefSeq, Apr 2018]

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To

be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u>. If you need a special design or shRNA sequence, please utilize our custom shRNA service.







## Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).