

## Product datasheet for **TR321484**

### CES3 Human shRNA Plasmid Kit (Locus ID 23491)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CES3 Human shRNA Plasmid Kit (Locus ID 23491)
Locus ID:	23491
Synonyms:	ES31
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	CES3 - Human, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 23491). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC053670</a> , <a href="#">NM_001185176</a> , <a href="#">NM_001185177</a> , <a href="#">NM_012122</a> , <a href="#">NM_024922</a> , <a href="#">NM_024922.3</a> , <a href="#">NM_024922.4</a> , <a href="#">NM_024922.5</a> , <a href="#">NM_001185177.1</a> , <a href="#">NM_001185176.1</a> , <a href="#">BC053670.1</a> , <a href="#">NM_024922.6</a> , <a href="#">NM_001185177.2</a>
UniProt ID:	<a href="#">Q6UWW8</a>
Summary:	This gene encodes a member of the carboxylesterase large family. The family members are responsible for the hydrolysis or transesterification of various xenobiotics, such as cocaine and heroin, and endogenous substrates with ester, thioester, or amide bonds. They may participate in fatty acyl and cholesterol ester metabolism, and may play a role in the blood-brain barrier system. This gene is expressed in several tissues, particularly in colon, trachea and in brain, and the protein participates in colon and neural drug metabolism. Multiple alternatively spliced transcript variants encoding distinct isoforms have been reported, but the biological validity and/or full-length nature of some variants have not been determined. [provided by RefSeq, Jun 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).