

Product datasheet for **TR500018**

Chrne Mouse shRNA Plasmid (Locus ID 11448)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Chrne Mouse shRNA Plasmid (Locus ID 11448)
Locus ID:	11448
Synonyms:	Ac; ACh; AChrepsilon; Acre; nAChRE
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Chrne - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 11448). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_009603 , NM_009603.1 , BC148476 , BC153108
UniProt ID:	P20782
Summary:	This gene encodes the epsilon subunit of the muscle-derived nicotinic acetylcholine receptor, a pentameric neurotransmitter receptor and member of the ligand-gated ion channel superfamily. The acetylcholine receptor changes subunit composition shortly after birth when the epsilon subunit replaces the gamma subunit seen in embryonic receptors. In mice, deficiency of this gene can lead to a decline in the number of nicotinic acetylcholine receptors at neuromuscular junctions and causes progressive muscle weakness, atrophy and premature death. Mutations in this gene serve as a pathophysiological model for human congenital myasthenia. Several alternatively spliced transcript variants of this gene have been described, but their full-length nature is not known. [provided by RefSeq, Nov 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).