

## Product datasheet for **TR500032**

### Ada Mouse shRNA Plasmid (Locus ID 11486)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Ada Mouse shRNA Plasmid (Locus ID 11486)
Locus ID:	11486
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Ada - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 11486). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC002075</a> , <a href="#">NM_001272052</a> , <a href="#">NM_007398</a> , <a href="#">NM_007398.1</a> , <a href="#">NM_007398.2</a> , <a href="#">NM_007398.3</a> , <a href="#">NM_007398.4</a> , <a href="#">NM_001272052.1</a>
UniProt ID:	<a href="#">P03958</a>
Summary:	Catalyzes the hydrolytic deamination of adenosine and 2-deoxyadenosine (PubMed:9272950). Plays an important role in purine metabolism and in adenosine homeostasis (PubMed:9272950, PubMed:10720488). Modulates signaling by extracellular adenosine, and so contributes indirectly to cellular signaling events (PubMed:11435465). Acts as a positive regulator of T-cell coactivation, by binding DPP4. Its interaction with DPP4 regulates lymphocyte-epithelial cell adhesion (By similarity). Enhances dendritic cell immunogenicity by affecting dendritic cell costimulatory molecule expression and cytokines and chemokines secretion (By similarity). Enhances CD4+ T-cell differentiation and proliferation (By similarity). Acts as a positive modulator of adenosine receptors ADORA1 and ADORA2A, by enhancing their ligand affinity via conformational change (By similarity). Stimulates plasminogen activation (By similarity). Plays a role in male fertility (By similarity). Plays a protective role in early postimplantation embryonic development (PubMed:9272950). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).