

## Product datasheet for **TR500052**

### **Adrb2 Mouse shRNA Plasmid (Locus ID 11555)**

#### **Product data:**

<b>Product Type:</b>	shRNA Plasmids
<b>Product Name:</b>	Adrb2 Mouse shRNA Plasmid (Locus ID 11555)
<b>Locus ID:</b>	11555
<b>Synonyms:</b>	Adrb; Adrb-2; Bad; Badm; Gpc; Gpcr7
<b>Vector:</b>	pRS (TR20003)
<b>E. coli Selection:</b>	Ampicillin
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Format:</b>	Retroviral plasmids
<b>Components:</b>	Adrb2 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 11555). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
<b>RefSeq:</b>	<a href="#">BC032883</a> , <a href="#">NM_007420</a> , <a href="#">NM_007420.1</a> , <a href="#">NM_007420.2</a> , <a href="#">NM_007420.3</a> , <a href="#">BC030346</a>
<b>UniProt ID:</b>	<a href="#">P18762</a>
<b>Summary:</b>	This intronless gene belongs to the G-protein-coupled receptor superfamily, which includes transmembrane proteins that play a role in signal transduction across biological membranes resulting in a variety of physiological responses. The encoded protein is a beta-2 adrenergic receptor which is activated by catecholamine ligands such as adrenaline and epinephrine. The protein participates in the classical signaling pathway involving G protein, adenylyl cyclase, cAMP (3'-5'-cyclic adenosine monophosphate) and protein kinase A (PKA). In humans, this gene is implicated in susceptibility to nocturnal asthma, obesity and type 2 diabetes. [provided by RefSeq, Apr 2013]
<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).