

## Product datasheet for **TR500231**

### Btrc Mouse shRNA Plasmid (Locus ID 12234)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Btrc Mouse shRNA Plasmid (Locus ID 12234)
Locus ID:	12234
Synonyms:	b-TrCP; beta-TrCP; Beta-Trcp1; E3RS-lkappaB; E3RSIkappaB; Fbw1a; FWD1; HOS; Slimb
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Btrc - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 12234). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<u><a href="#">BC003989</a></u> , <u><a href="#">NM_001037758</a></u> , <u><a href="#">NM_001286465</a></u> , <u><a href="#">NM_001286466</a></u> , <u><a href="#">NM_009771</a></u> , <u><a href="#">NM_001360120</a></u> , <u><a href="#">NM_001360124</a></u> , <u><a href="#">NM_009771.1</a></u> , <u><a href="#">NM_009771.2</a></u> , <u><a href="#">NM_009771.3</a></u> , <u><a href="#">NM_001037758.1</a></u> , <u><a href="#">NM_001037758.2</a></u> , <u><a href="#">NM_001286466.1</a></u> , <u><a href="#">NM_001286465.1</a></u> , <u><a href="#">BC052879</a></u> , <u><a href="#">NM_001360122</a></u> , <u><a href="#">NM_001360126</a></u> , <u><a href="#">NM_001360127</a></u>
UniProt ID:	<u><a href="#">Q3ULA2</a></u>



[View online »](#)

**Summary:** Substrate recognition component of a SCF (SKP1-CUL1-F-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins. Recognizes and binds to phosphorylated target proteins (PubMed:10097128, PubMed:16371461, PubMed:18782782, PubMed:9859996, PubMed:9990853, PubMed:21911472). SCF(BTRC) mediates the ubiquitination of phosphorylated NFKB, ATF4, CDC25A, DLG1, FBXO5, PER1, SMAD3, SMAD4, SNAI1 and probably NFKB2. SCF(BTRC) mediates the ubiquitination of CTNNB1 and participates in Wnt signaling (By similarity). SCF(BTRC) mediates the ubiquitination of NFKBIA, NFKBIB and NFKBIE; the degradation frees the associated NFKB1 to translocate into the nucleus and to activate transcription (PubMed:9859996, PubMed:10097128). Ubiquitination of NFKBIA occurs at 'Lys-21' and 'Lys-22' (PubMed:9859996, PubMed:10097128). SCF(BTRC) mediates the ubiquitination of CEP68; this is required for centriole separation during mitosis (By similarity). SCF(BTRC) mediates the ubiquitination and subsequent degradation of nuclear NFE2L1 (PubMed:21911472). Has an essential role in the control of the clock-dependent transcription via degradation of phosphorylated PER1 and PER2 (PubMed:18782782). May be involved in ubiquitination and subsequent proteasomal degradation through a DBB1-CUL4 E3 ubiquitin-protein ligase (By similarity). Required for activation of NFKB-mediated transcription by IL1B, MAP3K14, MAP3K1, IKBKB and TNF (By similarity). Required for proteolytic processing of GLI3 (PubMed:16371461). Mediates ubiquitination of REST, thereby leading to its proteasomal degradation (By similarity).[UniProtKB/Swiss-Prot Function]

**shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

**Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).