

## Product datasheet for **TR500512**

### Asap1 Mouse shRNA Plasmid (Locus ID 13196)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Asap1 Mouse shRNA Plasmid (Locus ID 13196)
Locus ID:	13196
Synonyms:	AV239055; Ddef1; DEF-1; mKIAA1249; PAP; s19
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Asap1 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector (Gene ID = 13196). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC094581</a> , <a href="#">NM_001276461</a> , <a href="#">NM_001276462</a> , <a href="#">NM_001276463</a> , <a href="#">NM_001276467</a> , <a href="#">NM_010026</a> , <a href="#">NM_010026.1</a> , <a href="#">NM_010026.2</a> , <a href="#">NM_001276461.1</a> , <a href="#">NM_001276467.1</a> , <a href="#">NM_001276463.1</a> , <a href="#">NM_001276462.1</a> , <a href="#">BC002201</a> , <a href="#">BC048818</a>
UniProt ID:	<a href="#">Q9QWY8</a>
Summary:	May function as a signal transduction protein involved in the differentiation of fibroblasts into adipocytes and possibly other cell types. Plays a role in ciliogenesis (By similarity). Possesses phosphatidylinositol 4,5-bisphosphate-dependent GTPase-activating protein activity for ARF1 (ADP ribosylation factor 1) and ARF5 and a lesser activity towards ARF6. May coordinate membrane trafficking with cell growth or actin cytoskeleton remodeling by binding to both SRC and PIP2.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).