

Product datasheet for **TR500539**

Dmpk Mouse shRNA Plasmid (Locus ID 13400)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Dmpk Mouse shRNA Plasmid (Locus ID 13400)
Locus ID:	13400
Synonyms:	DM; Dm15; DMK; MDPK; MT-PK
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Dmpk - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 13400). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC024150 , BC075715 , NM_001190490 , NM_001190491 , NM_032418 , NR_126095 , NM_001190491.1 , NM_032418.1 , NM_032418.2 , NM_001190490.1 , BC056615 , NM_001190490.2 , NM_001190491.2 , NM_032418.3
UniProt ID:	P54265
Summary:	The protein encoded by this gene is a serine/threonine protein kinase that contains coiled-coil and C-terminal membrane association domains. In the embryonic mouse, it is found in cardiac and skeletal myocytes where it appears to play a role in myogenesis. In adults, the transcript is localized to several tissues including brain, heart, and skeletal and smooth muscle, and a function in cytoskeletal remodeling has been described. Transcripts with expanded CUG repeats in the 3' untranslated region mediate alternative splicing of several genes and sequester RNA binding proteins and RNA transcripts that contain CAG repeats, resulting in myotonic dystrophy, an autosomal dominant neuromuscular disorder. Alternative splicing results in multiple protein coding and non-coding transcript variants. [provided by RefSeq, Oct 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).