

## Product datasheet for **TR500747**

### Gaa Mouse shRNA Plasmid (Locus ID 14387)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Gaa Mouse shRNA Plasmid (Locus ID 14387)
Locus ID:	14387
Synonyms:	E430018M07Rik
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Gaa - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 14387). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC010210</a> , <a href="#">NM_001159324</a> , <a href="#">NM_008064</a> , <a href="#">NM_001159324.1</a> , <a href="#">NM_001159324.2</a> , <a href="#">NM_008064.1</a> , <a href="#">NM_008064.2</a> , <a href="#">NM_008064.3</a> , <a href="#">NM_008064.4</a>
UniProt ID:	<a href="#">P70699</a>
Summary:	This gene encodes a lysosomal acid glucosidase that is involved in the degradation of glycogen. The encoded preproprotein undergoes proteolytic processing to generate a mature enzyme that cleaves alpha-1-4 and alpha-1-6 glycosidic bonds of glycogen, maltose and intermediate oligosaccharides within the lysosome. Mice lacking the encoded protein exhibit symptoms similar to human Pompe syndrome such as accumulation of glycogen in cardiac and skeletal muscle lysosomes resulting in reduced mobility and strength. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Nov 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).