

## Product datasheet for **TR501202**

### **Klra6 Mouse shRNA Plasmid (Locus ID 16637)**

#### **Product data:**

<b>Product Type:</b>	shRNA Plasmids
<b>Product Name:</b>	Klra6 Mouse shRNA Plasmid (Locus ID 16637)
<b>Locus ID:</b>	16637
<b>Synonyms:</b>	ly-49f; Ly49; Ly49f
<b>Vector:</b>	pRS (TR20003)
<b>E. coli Selection:</b>	Ampicillin
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Format:</b>	Retroviral plasmids
<b>Components:</b>	Klra6 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 16637). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
<b>RefSeq:</b>	<a href="#">NM_008464</a> , <a href="#">NM_008464.1</a> , <a href="#">NM_008464.2</a> , <a href="#">BC152859</a> , <a href="#">BC156612</a>
<b>UniProt ID:</b>	<a href="#">Q60653</a>
<b>Summary:</b>	This gene belongs to the highly polymorphic family of C-type lectin-like Ly49 genes that are expressed in natural killer (NK) cells. The encoded protein is a homodimeric type II transmembrane receptor located at the cell surface and inhibits NK cell activation upon ligand binding. This gene is located in a cluster of several Klra paralogs on chromosome 6. Different strains of mice show variation in the number of paralogs, including strain specific duplications, deletions and pseudogene sequences. [provided by RefSeq, Apr 2015]
<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).