

## Product datasheet for **TR501328**

### Mef2a Mouse shRNA Plasmid (Locus ID 17258)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Mef2a Mouse shRNA Plasmid (Locus ID 17258)
Locus ID:	17258
Synonyms:	A430079H05Rik
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Mef2a - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 17258). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC061128</a> , <a href="#">BC096598</a> , <a href="#">NM_001033713</a> , <a href="#">NM_001291191</a> , <a href="#">NM_001291192</a> , <a href="#">NM_001291195</a> , <a href="#">NM_001291196</a> , <a href="#">NM_001357324</a> , <a href="#">NM_001357325</a> , <a href="#">NM_001033713.1</a> , <a href="#">NM_001033713.2</a> , <a href="#">NM_001291196.1</a> , <a href="#">NM_001291195.1</a> , <a href="#">NM_001291191.1</a> , <a href="#">NM_001291192.1</a> , <a href="#">BC019116</a> , <a href="#">NM_001291191.2</a>
UniProt ID:	<a href="#">Q60929</a>
Summary:	Transcriptional activator which binds specifically to the MEF2 element, 5'-YTA[AT](4)TAR-3', found in numerous muscle-specific genes. Also involved in the activation of numerous growth factor- and stress-induced genes. Mediates cellular functions not only in skeletal and cardiac muscle development, but also in neuronal differentiation and survival. Plays diverse roles in the control of cell growth, survival and apoptosis via p38 MAPK signaling in muscle-specific and/or growth factor-related transcription. In cerebellar granule neurons, phosphorylated and sumoylated MEF2A represses transcription of NUR77 promoting synaptic differentiation. Associates with chromatin to the ZNF16 promoter (By similarity). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).