

## Product datasheet for **TR501333**

### Men1 Mouse shRNA Plasmid (Locus ID 17283)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Men1 Mouse shRNA Plasmid (Locus ID 17283)
Locus ID:	17283
Synonyms:	AW045611
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Men1 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 17283). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC036287</a> , <a href="#">NM_001168488</a> , <a href="#">NM_001168489</a> , <a href="#">NM_001168490</a> , <a href="#">NM_008583</a> , <a href="#">NM_001168490.1</a> , <a href="#">NM_001168488.1</a> , <a href="#">NM_008583.1</a> , <a href="#">NM_008583.2</a> , <a href="#">NM_001168489.1</a> , <a href="#">BC006610</a>
UniProt ID:	<a href="#">O88559</a>
Summary:	Essential component of a MLL/SET1 histone methyltransferase (HMT) complex, a complex that specifically methylates 'Lys-4' of histone H3 (H3K4). Functions as a transcriptional regulator. Binds to the TERT promoter and represses telomerase expression. Plays a role in TGFB1-mediated inhibition of cell-proliferation, possibly regulating SMAD3 transcriptional activity. Represses JUND-mediated transcriptional activation on AP1 sites, as well as that mediated by NFKB subunit RELA. Positively regulates HOXC8 and HOXC6 gene expression (By similarity). May be involved in normal hematopoiesis through the activation of HOXA9 expression. May be involved in DNA repair.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).