

## Product datasheet for **TR501994**

### Cxcl12 Mouse shRNA Plasmid (Locus ID 20315)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Cxcl12 Mouse shRNA Plasmid (Locus ID 20315)
Locus ID:	20315
Synonyms:	Pbsf; Scyb12; Sdf1; Tlsf; Tpar1
Vector:	pRS (TR20003)
Format:	Retroviral plasmids
Components:	Cxcl12 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 20315). 5µg purified plasmid DNA per construct Non-effective 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<u><a href="#">BC006640</a></u> , <u><a href="#">NM_001012477</a></u> , <u><a href="#">NM_013655</a></u> , <u><a href="#">NM_021704</a></u> , <u><a href="#">BC020297</a></u> , <u><a href="#">BC029335</a></u> , <u><a href="#">BC046827</a></u>
Summary:	This gene encodes a member of the alpha chemokine protein family. The encoded protein is secreted and functions as the ligand for the G-protein coupled receptor, chemokine (C-X-C motif) receptor 4. The encoded protein plays a role in many diverse cellular functions, including embryogenesis, immune surveillance, inflammation response, tissue homeostasis, and tumor growth and metastasis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).