

## Product datasheet for **TR502464**

### Zbtb7b Mouse shRNA Plasmid (Locus ID 22724)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Zbtb7b Mouse shRNA Plasmid (Locus ID 22724)
Locus ID:	22724
Synonyms:	c-Krox; Thpok; Zfp67
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Zbtb7b - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 22724). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC046533</a> , <a href="#">BC091651</a> , <a href="#">NM_009565</a> , <a href="#">NM_001355206</a> , <a href="#">NM_001355211</a> , <a href="#">NM_001355212</a> , <a href="#">NM_009565.1</a> , <a href="#">NM_009565.2</a> , <a href="#">NM_009565.3</a> , <a href="#">NM_009565.4</a> , <a href="#">BC089185</a> , <a href="#">NM_009565.5</a>
UniProt ID:	<a href="#">Q64321</a>



[View online »](#)

**Summary:**

Transcription regulator that acts as a key regulator of lineage commitment of immature T-cell precursors. Exerts distinct biological functions in the mammary epithelial cells and T cells in a tissue-specific manner (PubMed:15729333, PubMed:29420538). Necessary and sufficient for commitment of CD4 lineage, while its absence causes CD8 commitment. Development of immature T-cell precursors (thymocytes) to either the CD4 helper or CD8 killer T-cell lineages correlates precisely with their T-cell receptor specificity for major histocompatibility complex class II or class I molecules, respectively. Cross-antagonism between ZBTB7B and CBF complexes are determinative to CD4 versus CD8 cell fate decision (PubMed:15729333, PubMed:24880459, PubMed:18258917, PubMed:23481257). Suppresses RUNX3 expression and imposes CD4+ lineage fate by inducing the SOCS suppressors of cytokine signaling. induces, as a transcriptional activator, SOCS genes expression which represses RUNX3 expression and promotes the CD4+ lineage fate (PubMed:24880459). During CD4 lineage commitment, associates with multiple sites at the CD8 locus, acting as a negative regulator of the CD8 promoter and enhancers by epigenetic silencing through the recruitment of class II histone deacetylases, such as HDAC4 and HDAC5, to these loci (PubMed:22730529). Regulates the development of IL17-producing CD1d-restricted natural killer (NK) T cells (PubMed:23105140). Also functions as an important metabolic regulator in the lactating mammary glands. Critical feed-forward regulator of insulin signaling in mammary gland lactation, directly regulates expression of insulin receptor substrate-1 (IRS-1) and insulin-induced Akt-mTOR-SREBP signaling (PubMed:29420538). Transcriptional repressor of the collagen COL1A1 and COL1A2 genes. May also function as a repressor of fibronectin and possibly other extracellular matrix genes (PubMed:7937772). Potent driver of brown fat development, thermogenesis and cold-induced beige fat formation (PubMed:28784777). Recruits the brown fat lncRNA 1 (Blnc1):HNRNPU ribonucleoprotein complex to activate thermogenic gene expression in brown and beige adipocytes (PubMed:28784777). [UniProtKB/Swiss-Prot Function]

**shRNA Design:**

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).