

Product datasheet for **TR502678**

Trp53bp1 Mouse shRNA Plasmid (Locus ID 27223)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Trp53bp1 Mouse shRNA Plasmid (Locus ID 27223)
Locus ID:	27223
Synonyms:	53BP1; m53BP1; p53BP1; Tp53bp1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Trp53bp1 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 27223). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC035206 , BC079906 , NM_001290830 , NM_013735 , NM_013735.1 , NM_013735.2 , NM_013735.3 , NM_013735.4 , NM_001290830.1 , BC079906.1 , BC037378
UniProt ID:	P70399

[View online »](#)

Summary:	<p>Double-strand break (DSB) repair protein involved in response to DNA damage, telomere dynamics and class-switch recombination (CSR) during antibody genesis (PubMed:15159415, PubMed:15077110, PubMed:20453858, PubMed:23333305, PubMed:26308889, PubMed:20362325). Plays a key role in the repair of double-strand DNA breaks (DSBs) in response to DNA damage by promoting non-homologous end joining (NHEJ)-mediated repair of DSBs and specifically counteracting the function of the homologous recombination (HR) repair protein BRCA1 (PubMed:23333305, PubMed:20362325). In response to DSBs, phosphorylation by ATM promotes interaction with RIF1 and dissociation from NUDT16L1/TIRR, leading to recruitment to DSBs sites. Recruited to DSBs sites by recognizing and binding histone H2A monoubiquitinated at 'Lys-15' (H2AK15Ub) and histone H4 dimethylated at 'Lys-20' (H4K20me2), two histone marks that are present at DSBs sites. Required for immunoglobulin class-switch recombination (CSR) during antibody genesis, a process that involves the generation of DNA DSBs (PubMed:15159415, PubMed:15077110). Participates to the repair and the orientation of the broken DNA ends during CSR (PubMed:26308889). In contrast, it is not required for classic NHEJ and V(D)J recombination (PubMed:15159415). Promotes NHEJ of dysfunctional telomeres (By similarity). [UniProtKB/Swiss-Prot Function]</p>
shRNA Design:	<p>These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our custom shRNA service.</p>
Performance Guaranteed:	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>