

Product datasheet for **TR502705**

Abca7 Mouse shRNA Plasmid (Locus ID 27403)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Abca7 Mouse shRNA Plasmid (Locus ID 27403)
Locus ID:	27403
Synonyms:	Abc51; ABCX
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Abca7 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 27403). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001347081 , NM_013850 , NM_013850.1 , BC024511 , BC036566 , BC038521 , BC145120 , BC150708 , NM_013850.2
UniProt ID:	Q91V24
Summary:	The protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intracellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the ABC1 subfamily. Members of the ABC1 subfamily comprise the only major ABC subfamily found exclusively in multicellular eukaryotes. This protein is widely expressed with highest detection in spleen and hematopoietic tissues. Defects in this gene cause an increase in amyloid-beta deposits in a mouse model of Alzheimer's disease, and a related human protein is thought to play a role in lipid homeostasis in cells of the immune system. [provided by RefSeq, Jan 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).