

## Product datasheet for **TR502831**

### Mapk6 Mouse shRNA Plasmid (Locus ID 50772)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Mapk6 Mouse shRNA Plasmid (Locus ID 50772)
Locus ID:	50772
Synonyms:	2610021I23Rik; D130053K17Rik; Erk3; Mapk4; Mapk63; Prkm4; Prkm6
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Mapk6 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 50772). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC048779</a> , <a href="#">BC052420</a> , <a href="#">BC100385</a> , <a href="#">NM_015806</a> , <a href="#">NM_027418</a> , <a href="#">NM_015806.1</a> , <a href="#">NM_015806.2</a> , <a href="#">NM_015806.3</a> , <a href="#">NM_015806.4</a> , <a href="#">NM_015806.5</a> , <a href="#">NM_027418.1</a> , <a href="#">NM_027418.2</a> , <a href="#">BC006778</a> , <a href="#">BC024684</a> , <a href="#">BC037060</a>
UniProt ID:	<a href="#">Q61532</a>
Summary:	Atypical MAPK protein. Phosphorylates microtubule-associated protein 2 (MAP2) and MAPKAPK5. The precise role of the complex formed with MAPKAPK5 is still unclear, but the complex follows a complex set of phosphorylation events: upon interaction with atypical MAPKAPK5, ERK3/MAPK6 is phosphorylated at Ser-189 and then mediates phosphorylation and activation of MAPKAPK5, which in turn phosphorylates ERK3/MAPK6. May promote entry in the cell cycle.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).