

## Product datasheet for **TR502934**

### Ube2k Mouse shRNA Plasmid (Locus ID 53323)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Ube2k Mouse shRNA Plasmid (Locus ID 53323)
Locus ID:	53323
Synonyms:	AW492011; D5Ertd601e; E2-25k; HIP-2; Hip2; Hypg; Lig
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Ube2k - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector (Gene ID = 53323). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC002013</a> , <a href="#">BC085311</a> , <a href="#">NM_001310618</a> , <a href="#">NM_001310619</a> , <a href="#">NM_016786</a> , <a href="#">NR_144566</a> , <a href="#">NM_016786.1</a> , <a href="#">NM_016786.2</a> , <a href="#">NM_016786.3</a> , <a href="#">NM_016786.4</a>
UniProt ID:	<a href="#">P61087</a>
Summary:	Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In vitro, in the presence or in the absence of BRCA1-BARD1 E3 ubiquitin-protein ligase complex, catalyzes the synthesis of 'Lys-48'-linked polyubiquitin chains. Does not transfer ubiquitin directly to but elongates monoubiquitinated substrate protein. Mediates the selective degradation of short-lived and abnormal proteins, such as the endoplasmic reticulum-associated degradation (ERAD) of misfolded luminal proteins. Ubiquitinates huntingtin. May mediate foam cell formation by the suppression of apoptosis of lipid-bearing macrophages through ubiquitination and subsequent degradation of p53/TP53. Proposed to be involved in ubiquitination and proteolytic processing of NF-kappa-B; in vitro supports ubiquitination of NFKB1. Involved in stabilization of CASP12 during ER stress-mediated amyloid-beta neurotoxicity probably by inhibiting proteasome activity; in vitro ubiquitinates CASP12.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).