

Product datasheet for **TR503011**

Atf7ip Mouse shRNA Plasmid (Locus ID 54343)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Atf7ip Mouse shRNA Plasmid (Locus ID 54343)
Locus ID:	54343
Synonyms:	2610204M12Rik; AM; Mcaf1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Atf7ip - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 54343). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC052460 , BC060246 , NM_019426 , NM_019426.1 , NM_019426.2 , BC022714 , BC037617 , BC048827
UniProt ID:	Q7TT18
Summary:	Recruiter that couples transcriptional factors to general transcription apparatus and thereby modulates transcription regulation and chromatin formation. Can both act as an activator or a repressor depending on the context (PubMed:10777215). Required for HUSH-mediated heterochromatin formation and gene silencing (By similarity). Mediates MBD1-dependent transcriptional repression, probably by recruiting complexes containing SETDB1. Stabilizes SETDB1, is required to stimulate histone methyltransferase activity of SETDB1 and facilitates the conversion of dimethylated to trimethylated H3 'Lys-9' (H3K9me3). The complex formed with MBD1 and SETDB1 represses transcription and couples DNA methylation and histone H3 'Lys-9' trimethylation (H3K9me3) (PubMed:14536086). Facilitates telomerase TERT and TERC gene expression by SP1 in cancer cells (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).