

Product datasheet for **TR503155**

Arl3 Mouse shRNA Plasmid (Locus ID 56350)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Arl3 Mouse shRNA Plasmid (Locus ID 56350)
Locus ID:	56350
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Arl3 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 56350). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC042941 , NM_019718 , NM_001355233 , NM_019718.1 , NM_019718.2 , NM_019718.3
UniProt ID:	Q9WUL7
Summary:	Small GTP-binding protein which cycles between an inactive GDP-bound and an active GTP-bound form, and the rate of cycling is regulated by guanine nucleotide exchange factors (GEF) and GTPase-activating proteins (GAP) (PubMed:18376416). Required for normal cytokinesis and cilia signaling. Required for targeting proteins to the cilium, including myristoylated NPHP3 and prenylated INPP5E. Targets NPHP3 to the ciliary membrane by releasing myristoylated NPHP3 from UNC119B cargo adapter into the cilium (By similarity). Requires assistance from GTPase-activating proteins (GAPs) like RP2 and PDE6D, in order to cycle between inactive GDP-bound and active GTP-bound forms (PubMed:15979089). Required for PKD1:PKD2 complex targeting from the trans-Golgi network to the cilium (PubMed:25405894).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).