

## Product datasheet for **TR503189**

### Hbs1l Mouse shRNA Plasmid (Locus ID 56422)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Hbs1l Mouse shRNA Plasmid (Locus ID 56422)
Locus ID:	56422
Synonyms:	2810035F15Rik; AI326327; eRFS
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Hbs1l - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 56422). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">BC010251</a> , <a href="#">BC085235</a> , <a href="#">NM_001042593</a> , <a href="#">NM_001145209</a> , <a href="#">NM_019702</a> , <a href="#">NM_019702.1</a> , <a href="#">NM_019702.2</a> , <a href="#">NM_001042593.1</a> , <a href="#">NM_001145209.1</a> , <a href="#">NM_001145209.2</a>
UniProt ID:	<a href="#">Q69ZS7</a>
Summary:	Cotranslational quality control factor involved in the No-Go Decay (NGD) pathway. In the presence of ABCE1 and PELO, is required for 48S complex formation from 80S ribosomes and dissociation of vacant 80S ribosomes. Together with PELO and in presence of ABCE1, recognizes stalled ribosomes and promotes dissociation of elongation complexes assembled on non-stop mRNAs; this triggers endonucleolytic cleavage of the mRNA, a mechanism to release non-functional ribosomes and to degrade damaged mRNAs as part of the No-Go Decay (NGD) pathway.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).