

## Product datasheet for **TR503398**

### Huwe1 Mouse shRNA Plasmid (Locus ID 59026)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Huwe1 Mouse shRNA Plasmid (Locus ID 59026)
Locus ID:	59026
Synonyms:	5430439H10Rik; Arf-bp1; AU041296; C80292; C430014N20Rik; Gm1718; Ib772; LASU1; Mule; Ureb1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Huwe1 - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 59026). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_021523</a> , <a href="#">BC004007</a> , <a href="#">BC011391</a> , <a href="#">BC017642</a> , <a href="#">BC054372</a> , <a href="#">BC060979</a> , <a href="#">BC070444</a> , <a href="#">BC079665</a>
Summary:	E3 ubiquitin-protein ligase which mediates ubiquitination and subsequent proteasomal degradation of target proteins. Regulates apoptosis by catalyzing the polyubiquitination and degradation of MCL1. Mediates monoubiquitination of DNA polymerase beta (POLB) at 'Lys-41', 'Lys-61' and 'Lys-81', thereby playing a role in base-excision repair. Also ubiquitinates the p53/TP53 tumor suppressor and core histones including H1, H2A, H2B, H3 and H4. Binds to an upstream initiator-like sequence in the preprodynorphin gene (By similarity). Regulates neural differentiation and proliferation by catalyzing the polyubiquitination and degradation of MYCN. May regulate abundance of CDC6 after DNA damage by polyubiquitinating and targeting CDC6 to degradation. Mediates polyubiquitination of PA2G4 (By similarity). Acts in concert with MYCBP2 to regulate the circadian clock gene expression by promoting the lithium-induced ubiquitination and degradation of NR1D1 (By similarity).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).