

Product datasheet for **TR503638**

Pyurf Mouse shRNA Plasmid (Locus ID 66459)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Pyurf Mouse shRNA Plasmid (Locus ID 66459)
Locus ID:	66459
Synonyms:	2610022G08Rik; AI847956; Prey
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Pyurf - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector (Gene ID = 66459). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	BC029232 , NM_025574 , NM_025574.1 , NM_025574.2 , NM_025574.3
UniProt ID:	Q9D1C3
Summary:	This gene encodes a small protein with a conserved DUF343 domain. The human ortholog of this gene expresses two distinct proteins from upstream and downstream coding regions. The upstream CDS encoding a DUF343 domain-containing protein has been conserved at this mouse locus, but the downstream CDS encoding a subunit of an enzyme involved in glycosylphosphatidylinositol biosynthesis has not been conserved. Instead, a separate locus on mouse chromosome 9 encodes the mouse homolog of the human phosphatidylinositol glycan anchor biosynthesis, class Y protein. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).